MANUALE D'USO
USER MANUAL
MANUEL D'UTILISATION
BENUTZERHANDBUCH
GEBRUIKERSHANDLEIDING

pag.3

pag.39

pag.75

s.111

pag.147

ARIA 1.0 Mg

ARIA 2.0 Al

ARIA ULTRA

ARIA KID

ARIA Speciale

NRIN®

www.ariawheels.com

INDICE

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INDEX

1. INTRODUCTION	pg.40
1.1 Information related to this user's manual	
1.2 Explanation of symbols	
1.3 After Sale warranty	
1.4 Use in compliance with the expected purposes	
1.5 Limitation of liability	
1.6 Duration of the wheelchair	
2. SAFETY	pg.42
2.1 Safety information	
2.2 Contraindications	
2.3 Safety equipment	
3. STRUCTURE AND OPERATION	pg.43
3.1 Components of the wheelchair and materials	
3.2 Parking brakes	
3.3 Backrests	
3.4 Push handles	
3.5 Mudguard	
3.6 Footrest	
3.7 Anti-tip devices	
3.8 Centre of gravity	
3.9 Padded seat cushion	
3.10 Camber	
3.11 Sitting balance	
4 ACCESSORIES	pg.51
5 PUT INTO SERVICE	pg.52
5.1 Safety indications	
6 MOVING WITH THE WHEELCHAIR	pg.52
6.1 Safety indications	
6.2 Braking during movement	
6.3 Sitting down and getting up from the wheelchair	
6.4 Forward movement and changing direction with the wheelchair	
6.5 Going up and down steps and drops	
6.6 Going up and down stairs	
6.7 Going up ramps and sloping surfaces	
6.8 Stability and balance	
7. TRANSPORTATION	pg.60
7.1 Safety indications	
7.2 Closing and opening the wheelchair	
7.3 Assembling and dismantling of wheels	
7.4 Transportation of the occupied wheelchair in a vehicle	
8. MAINTENANCE	pg.66
8.1 Safety indications	
8.2 Maintenance plan	
8.3 Care	
8.4 Disinfection	
9. BREAKDOWN SOLUTIONS	pg.70
9.1 Safety indications	
9.2 Identifying and solving breakdowns	
10. AFTER USE	pg.72
10.1 Reuse	
10.2 Disposal	
11. TECHNICAL DATA	pg.72
11.1 Weight and dimensions	
11.2 Environmental conditions	
11.3 Tyres. The optimal pressure depends on the type of tyre	
11.4 Labels	
12. DECLARATION OF CONFORMITY	pg.74

1 INTRODUCTION

1.1 INFORMATION RELATED TO THE PRESENT USER'S MANUAL

We thank you for choosing the ARIA wheelchair.

This user's manual contains important indications for the use of the wheelchair. To be able to utilise the wheelchair and move safely, we ask you to read the manual carefully and to observe the safety indications.

In the event of vision problems, it is possible to view this user's manual as a PDF file on the Internet at www.ariawheels.com and zoom in on the screen according to your needs. If you are unable to enlarge the texts and graphics enough, please contact an ARIA official distributor in your national territory; the addresses are listed on the website. If necessary, we will provide a high resolution PDF file of the manual. Additionally, it is possible to obtain the reading of the PDF file using appropriate software that use special linguistic functions installed on your computer.

1.2 EXPLANATION OF SYMBOLS

Symbols in this manual

All the instructions in this manual are labelled with symbols. In front of the individual messages are symbols and words that indicate the gravity of the danger.

WARNING !

Indicates a potential danger that could cause severe injuries or death, if not avoided.

WARNING !

CAUTION !

Indicates a potential danger that could cause small to medium injuries, if not avoided.

CAUTION !

IMPORTANT !

Indicates a potential danger that could cause material damages, if not avoided. Indicates tips and useful recommendations for an efficient and inconvenience-free use of the wheelchair.

IMPORTANT !

This product is in compliance with the directive (93/42/ECC) regarding medical products. The launch date of this product is indicated in the declaration of compliance CE.

Symbols on the product

An identification label is applied on the product.

The symbols on the label are defined in detail in section 11.4 of this user's manual.

1.3 AFTER SALE WARRANTY

Aria Wheels Srl guarantees the absence of defects and the functionality of its wheelchairs. The warranty covers defects that arise from manufacturing, the use of bad quality materials or the inappropriate working of these. The warranty rights towards Aria Wheels Srl can only be executed by the seller, and not the user of the wheelchair. The commercial warranty does not cover the normal wear, damages or faults that come directly or indirectly from accidents, falls, collisions, improper use, insufficient maintenance and alterations made by non authorised personnel.

The parts subject to wear and particular consumption (tyres, seat canvases and backrest...) are not covered by warranty. The warranty falls if any unexpected alterations are made to the wheelchair, or if alterations are made with inadequate or unoriginal parts. In case of not authorized alterations, the warranty does not cover the costs arising from the elimination of the defect such as shipping costs, labour costs, expenses etc.

The warranty has a duration of 24 (twenty-four) months starting from the date of invoice.

All the components not manufactured by Aria Wheels Srl have a warranty from their own manufacturer.

Moreover, the terms and conditions are an integral part of the general and specific conditions for the individual countries in which the product is commercialised.

1.4 USE IN COMPLIANCE WITH THE EXPECTED PURPOSES.

The ARIA wheelchair is to be pushed manually and serves the sole purpose of movement of a person affected by motor disabilities, either autonomously or with a helper. Without a helper it should only be utilised by people who are physically and mentally capable to control and drive the wheelchair safely (for example moving forwards, changing direction and braking). The wheelchair should be used exclusively on flat surfaces, accessible turf and indoors. This active wheelchair has to be approved and appropriate for the specific physical conditions of the user.

WARNING!

Every use not compliant to the expected purposes or of any other kind can lead to dangerous situations.

1.5 LIMITATION OF LIABILITY

Aria Wheels Srl accepts no liability for damage arising from:

Non-compliance with the User Manual; Incorrect use; Natural wear and tear; Incorrect assembly or set-up by the purchaser or a third party; Technical modifications; Unauthorised modifications and/or use of unsuitable spare parts.

A written authorisation of Aria Wheels Srl must be obtained before installing additional adaptations on a ARIA wheelchair. Otherwise no liability claims can be made.

1.6 DURATION OF THE WHEELCHAIR

The expected duration of the product is four years, assuming a daily use at the condition that is used in compliance with the safety, maintenance and expected use instructions indicated in this manual.

2. SAFETY

2.1 SAFETY INFORMATION

This chapter contains important safety aspects, aimed at protecting the user of the wheelchair and his helper (when present), as well as guaranteeing a safe and inconvenience-free use of the wheelchair.

WARNING!

*Danger of accidents and severe injuries to the patient
If the wheelchair is not well adjusted it can lead to accidents with severe injuries.
The adjustments to the wheelchair should always be made by a specialised dealer.*

WARNING!

*Danger arising from inadequate driving behaviour
Moving on wet surfaces, gravel and uneven ground can lead to slipping.
Always adapt the speed and driving behaviour to the different situations
(atmospheric conditions, surface, individual abilities, etc.).*

WARNING!

*Danger of injuries to patient. In the event of a collision, it is possible to be wounded in areas that stick out of the wheelchair (such as feet or hands).
Avoid all collisions without having tried to use the brakes.
Never drive the wheelchair into an object.
Drive through narrow passages with caution.*

WARNING!

*Danger due to profligate driving.
Moving at high speeds you can lose control of your wheelchair and fall.
Never exceed the speed of 5km/h.
Avoid all types of collision.*

CAUTION!

*Danger of burning.
The components of the wheelchair can heat up if exposed to strong solar radiations.*

2.2 CONTRAINDICATIONS

No contraindications were found except in cases of proven sensitivity to the component materials.

CAUTION!

Risk of squashing your fingers!

There is always the risk of squashing your fingers or arms in moving parts of the wheelchair. Activating the folding mechanisms or inserting moving parts, such as the extractable shaft of the rear wheel, or the anti-tip device, make sure that nothing can get stuck.

2.3 SAFETY EQUIPMENT

WARNING!


Danger of accident!

Safety equipment (brakes, anti-tip device), if wrongly adjusted or not working, can cause accidents. Before each use of the wheelchair check the functioning of the safety equipment and have them regularly checked by a specialised dealer. specializzato.

The operation of safety equipment is described on chapter 3 Structure and operations.

3 STRUCTURE AND OPERATIONS

3.1 Components of the wheelchair and materials

ITEM	MATERIAL
Postural backrest and cover	Backrest body: Carbon Fiber or Aluminum Cover: breathable fabric and polyurethane foam
Mudguard	Carbon Fiber or ABS
Armrest	Armrest structure: Aluminum Integrated mudguard: Carbon Fiber Padding: Polyurethane (PU)
Frame	Aria 2.0Al: Aluminium alloy Aria 1.0Mg: Magnesium alloy ULTRA: Magnesium alloy SPECIALE: Magnesium alloy
Footrest	Aluminium alloy or Carbon Fiber
Tension back	Strap supports: Polyurethane and Fabric Band: Fabric and Velcro Cover: breathable fabric, polyurethane foam and leather
Seat	Structure: Carbon Fiber or Upholstery Cushion: not supplied

AVAILABLE WHEELCHAIR MODELS

Brief description

The Aria Wheels Srl wheelchair models are of the superlight kind with a rigid frame.

"The effective equipment of the wheelchair can vary from the pictures in the manual, as each wheelchair is built according to the indications given in the order."

3.2 PARKING BRAKES

The parking brakes are made to stop the wheelchair and impede any involuntary movements.
We use Push to Lock or Scissor Lock types.

WARNING!

Risk of falling following a sudden stop!

If during movement you use the parking brakes, it is no longer possible to control the direction and the wheelchair could become blocked and cause collisions or a fall from the wheelchair. Never pull the parking brakes whilst moving.



WARNING!

Danger of falling!



The functioning of the parking brakes is only guaranteed if the tyres have sufficient pressure.

Check the correct pressure of the tyres, 11.3 Tyres.


PUSH TO LOCK BRAKE

SCISSOR LOCK BRAKE

ULTRA'S BRAKE

3.3 BACKREST

Backrest cover

Aria wheelchair is supplied with a padded cover for the backrest in carbon fibre.
The cover in breathable material is removable thanks to an ad-hoc zipper and conveniently washable


Height of the postural backrest

It is possible to adjust the height of the backrest. Remove the upholstery, unscrew the 4 screws and screw them back in at the desired height. The adjustment has to be made by the specialised dealer.


ULTRA / SPECIALE / KID tension backrest

The tension of the backrest straps can be adjusted. Remove the upholstery, tension the bands according to the desired posture. The adjustment has to be made by the specialised dealer.

1.0 Mg


ULTRA 2.0 Al


Inclination of the postural backrest

It is possible to adjust the inclination of the backrest unscrewing the two screws on the clamp, adjusting the inclination and screwing them back in. The adjustment should be done by the specialised dealer.


3.4 PUSH HANDLES

The ARIA 1.0 Mg and Air 2.0 Al wheelchairs can be equipped with push bar (that also work as a bag holder) that are placed inside the clamp for adjusting the inclination of the backrest. The ARIA ULTRA wheelchair can be equipped with push handles. The height of both the push bar and the push handles can be adjusted. Unscrew the screws and screw them back in at the desired height. The adjustment has to be made by the specialised dealer.


1.0  2.0 


ULTRA

3.5 MUDGUARD

Depending on the request, the mudguard can be straight or with the tire cover curve. The straight model is also mounted in combination with the armrests. Both models are exchangeable with the exception that the curve model cannot be mounted in combination with the armrests.

1.0  2.0 


ULTRA


3.6 FOOTREST

The height of the footrest can be adjusted by unscrewing the two screws on the frame, adjusting the desired height and tightening it again. Also it is possible to adjust the angle of the footrest by loosening the screws under the aluminium or carbon fibre plate. Adjusted the angle to retighten the previously loosened screws. Adjustments must be performed by the specialist dealer.

1.0 Mg


2.0 Al

ULTRN

3.7 ANTI-TIP DEVICE

An anti-tip device stops the wheelchair from tipping backwards. The device is adjustable based on different needs. It is fixed to the central element in carbon fibre thanks to a clamp. It can be positioned on the right or left side

WARNING!

Danger of fall!

A wrongly adjusted or non-working anti-tip device can cause falls.

Before every use of the wheelchair control it is functioning and, if necessary, have the anti-tip device fitted or adjusted by a specialised dealer.

WARNING!

Danger of tipping

On uneven or soft ground, the anti-tip device can sink in holes or in the ground, thus limiting or compromising its correct safety function.


Only use the anti-tip device when you are on a flat and compact surface.

3.8 CENTER OF GRAVITY

Except for the model SPECIALE, It is possible to adjust the center of gravity unscrewing the four screws on the framework (two per side), adjusting the desired center of gravity and screwing them back in. The adjustment should be done by the specialised dealer.

1.0 Mg **2.0** Al

ULTRON

3.9 PADDED SEAT CUSHION

An appropriate padded cushion is necessary to guarantee an optimal distribution of pressure on the seat. The cushion is fixed to the seat stopping it from slipping thanks to Velcro.

The padded cushion is considered an option of the base model.

3.10 CAMBER

The wheels camber is selected when ordering the ARIA wheelchair. Subsequently it can be modified by replacing the wheels axle. The replacing should be made by the specialised dealer.

3.11 SITTING BALANCE

Anti-tip stability

The seating position, therefore the anti-tip stability, can be adjusted. The adjustment should be made by the specialised dealer.

Front height of seat


The front height of the seat can be adjusted. The adjustment should be made by the specialised dealer.

Back height of seating


The back height of the seat can be adjusted of a few centimetres.

The adjustment should be made by the specialised dealer.

1.0  2.0 

ULTRON


4. ACCESSORIES

ANTI-TIP AND TIPPING LEVER

The tipping lever makes it easy for an assistant to tilt the wheelchair backwards in order to go over a step. To do this turn the anti tipper down using your foot (as in the picture on the right).

PUSH HANDLE (AND BAG HOLDER)

5. PUT INTO SERVICE

5.1 SAFETY INDICATIONS

CAUTION!

*Danger of injuries to patient.
Before putting the wheelchair into service check it is in a good state and the functioning of the most important functions, See the maintenance plan, sec. 8.*

Your trusted dealer will prepare the wheelchair for use. The dealer will explain the main functions and will make sure that the wheelchair matches your requests and needs.

6. MOVING WITH THE WHEELCHAIR

6.1 SAFETY INDICATIONS

WARNING!

*Danger of accident!
The movement behaviour of the wheelchair could be considerably compromised if the pressure of the tyres is not the same.
Check the pressure of the tyres before any movement.*

CAUTION!

*Danger of crushing!
The distance between the rear wheel and the wheel's mudguard could be so small to cause the crushing of fingers.
Always make sure to move the wheelchair using solely the handrims.*

CAUTION!

*Danger of crushing!
The distance between the rear wheel and the parking brake could be so small to cause the crushing of fingers.
Always make sure to move the wheelchair using solely the handrims.*

6.2 BRAKING DURING MOVEMENT

During movement it is possible to brake the wheelchair putting pressure with your hands on the handrims.

WARNING!

Risk of falling after sudden braking!

If you use the parking brakes during movement, it is no longer possible to control the direction of the movement and the wheelchair could become suddenly blocked, which could lead to a collision or fall from the wheelchair. Never pull the parking brakes during movement.

WARNING!

Risk of falling off the wheelchair If the wheelchair is quickly decelerated by a helper who pulls the push handles, the user could fall from the wheelchair. Always fasten the seatbelt if present.

Make sure that the helper received individual training on how to push occupied wheelchairs.

CAUTION!

Risk of burns on the skin

The handrims could warm up, if you brake the wheelchair for a long period of time.

You could burn the palm of your hands.

Wear an appropriate pair of gloves.

Grab the handrims and exercise a homogenous pressure with both hands, until the wheelchair stops completely.

6.3 SITTING DOWN AND GETTING UP FROM THE WHEELCHAIR

WARNING!

Risk of falling!

During transfers the risk of falling is very high.

Only sit down and get up from the wheelchair with no help if you are physically able to do so.

WARNING!

Risk of falling!

*The wheelchair could tip forwards if you are standing on the footrest.
Never step on the footrest whilst sitting down or getting up from the wheelchair.*

CAUTION!

*If the brakes are disarmed or broken, the wheelchair could move in an uncontrolled manner.
Do not lean on the brakes whilst sitting down or standing up from the wheelchair.*

1. Pull the parking brakes.
2. Place your feet on the ground.
3. Hold on tightly to the wheelchair and, if necessary, to a fixed object nearby.
4. Slowly move onto the chair.

6.4 FORWARD MOVEMENT AND CHANGING DIRECTION WITH THE WHEELCHAIR

Using the handrims it is possible to move the wheelchair forwards and change direction.
Before moving without a helper, it is necessary to identify the tipping point of the wheelchair.

WARNING!

Danger of tipping

*The wheelchair could tip backwards, if an anti-tip device is not installed.
When determining the tipping point of the wheelchair, a helper should be placed directly behind the wheelchair to catch the wheelchair before it tips.
To stop the wheelchair from tipping, install an anti-tip device.*

WARNING!

Danger of tipping

*The wheelchair could tip forwards.
When adjusting the wheelchair, check its behaviour in the event of forward tipping and adjust your driving style accordingly.*

CAUTION!

Danger of tipping

A heavy load attached to the backrest could alter the centre of gravity of the wheelchair.

In this event, adjust your driving behaviour to account for the load.

Determining the tipping point

1. Release the brake.
2. Briefly move forwards, hold onto both handrims and move backwards giving a slight push.
3. Determine the tipping point by shifting your weight and contrasting the movement of the handrims. 2

6.5 GOING UP AND DOWN STEPS AND DROPS

WARNING!

Danger of falling!

Going up and down steps you could lose balance and tip the wheelchair.

Always move slowly and carefully when going over steps, such as pavements, and drops. Do not go over steps taller than 25cm.

CAUTION!

An anti-tip device prevents the wheelchair from tipping backwards.

Take off the anti-tip device before going over steps or drops.

Going down a step with a helper

1. Bring the wheelchair to the edge of the step and grab the push handle
2. The helper firmly grabs the push handles, places a foot on the anti-tip device (if installed) and inclines the wheelchair backwards, so that the front wheels lift from the ground.
3. The helper keeps the wheelchair in this position, pushes it carefully over the step and inclines the wheelchair forwards so that the front wheels are touch the ground again.

Going down a step without a helper

WARNING!

Danger of tipping

If you go down a step without a helper, there is a risk a backward tipping, if you do not have control of the wheelchair.

Initially learn how to go down a step with the aid of a helper.

Learn how to balance on the back wheels.

1. Bring the wheelchair to the edge of the step, lift the front wheels and balance.
2. Let the back wheels slowly slide along the edge. Firmly grasp with both hands the handrims and hold on until the front wheels are touching the ground again.


WARNING!

Risk of injuries to the user.

The backrest of the wheelchair can break earlier than expected if you frequently go up steps and pavements, it is suggested take the wheelchair from the frame or from the push handle if installed. The user could fall from the wheelchair. Take off the anti-tip device when going up steps and pavements.

1. Push the wheelchair backwards with the back wheels until the edge of the step.
2. The helper grabs the push handles and inclines the wheelchair so that the front wheels are off the ground; he then pulls the back wheels over the edge of the step until it is possible to put the front wheels back on the ground.

6.6 GOING UP AND DOWN STAIRS

WARNING!

Danger of falling.

Going up and down stairs you could lose balance and fall off the wheelchair. Always go up and down stairs with more than one step only with 2 helpers.

1. It is possible to go down flights of stairs one step at a time, as described in the picture. The first helper will be behind the wheelchair and will hold onto the push handles. The second helper will hold onto a fixed part of the front framework, ensuring the position of the wheelchair is correct from the front.

6.7 GOING UP RAMPS AND SLOPING SURFACES

WARNING!

Danger due to an uncontrolled movement! Going through sloping surfaces, the wheelchair could tip backwards, sideways or forwards. Go through long sloping surfaces essentially with a helper behind the wheelchair. Avoid leaning sideways. Avoid surfaces with a gradient bigger than 7°. Avoid abruptly changing direction on a sloping surface.


CAUTION!

The wheelchair could involuntarily move even on slightly sloping surfaces, if not controlled using the handrims. Apply the parking brakes if on a sloping surface with the wheelchair.

Uphill

To move uphill you need to give a slight push, keep up the momentum and at the same time control the direction.


- 1. Bend forwards and move the wheelchair forwards giving quick and energetic pushes to both handrims.*

Downhill

Moving downhill, it is important to control the direction of movement and, above all, the speed.

1. Lean against the backrest and carefully let the handrims slide through your hands. You must always be able to stop the wheelchair using the handrims.

WARNING!

Danger of burns

The handrims can warm up, if you brake the wheelchair for a long period of time. You could burn the palm of your hands.

Wear an appropriate pair of gloves.

6.8 STABILITY AND BALANCE

Some activities and movements during daily use can require you to stick out of the wheelchair, leaning forwards, sideways or backwards. This considerably affects the stability of the wheelchair. To be sure to always maintain balance, bear in mind the following:

Leaning forwards

WARNING!

Danger of falling off the wheelchair!

If you lean forwards you could fall off the wheelchair.

Do not bend forwards too much and do not move your body beyond the seat to reach for an object.

Do not bend forwards widening your knees to reach for an object on the ground.

1. Align the front wheels (to do so move the wheelchair slightly forwards and then backwards).
2. Apply the two parking brakes.
3. Lean forwards only up to the point where the upper body is still above the front wheels.

Leaning backwards

WARNING!

Danger of falling from the wheelchair!

If you lean backwards too far, you could tip the wheelchair.

Don't lean backwards beyond the wheelchair.

Use an anti-tip device.

1. Align the front wheels (to do so move the wheelchair slightly forwards and then backwards)
2. Do not apply the parking brakes.
3. Only grab the reachable objects without having to modify your sitting position.

7. TRANSPORTATION

7.1 SAFETY INDICATIONS

WARNING!

Risk of injuries if the wheelchair is not correctly and securely fixed

In the event of an accident, brake, etc. there can be severe injuries due to parts of the wheelchair not fixed. Always extract the wheels when transporting the wheelchair.

Firmly secure all the components of the wheelchair in the vehicle in order to stop them from loosening during transportation.



CAUTION!

Excessive abrasion of the material could compromise the resistance of the supporting parts. Don't move the wheelchair with the wheels dismantled on an abrasive surface (e.g. pulling the framework on tarmac).

7.2 CLOSING AND OPENING THE WHEELCHAIR

The backrest is to be reclined pulling the cord in the rear part immediately under the seat, then pull the backrest forwards and lower it onto the seat. (Only for Aria 1.0 Mg and Aria 2.0 Al)


WARNING!

Danger of tipping!

If the quick-release axle of the rear wheel is not completely blocked, the wheel can fall off during use. This can cause tipping. Always verify that the quick-release axles are completely blocked every time you assemble a wheel.


Dismantling of rear wheels

Press on the quick-release pin and extract the wheel at the same time.


Removing the casters

Unscrew the axle shaft with a special Allen key and remove the wheel.


Assembling of rear wheels

Press on the quick-release pin and position the wheel at the same time until the pin blocks it.

Assembly the casters

Follow the disassembly process on the contrary, taking care to tighten the axle pin securely.

7.4 TRANSPORTATION OF THE OCCUPIED WHEELCHAIR IN A VEHICLE.

Even when the ARIA wheelchair is correctly fixed and the following rules are followed, there could be injuries to the passengers in the event of an impact or an abrupt stop. Therefore, Aria Wheels Srl strongly recommends to transfer the user of the wheelchair onto a seat of the vehicle and to fasten his seatbelt.

Do not make any modifications or substitutions to the wheelchair (structure, framework or parts) without a written consent by Aria Wheels Srl

WARNING!

Risk of sever injuries or death

To use the wheelchair as a seat in the vehicle the height of the backrest must be of at least 400mm.

For the transportation of the occupied wheelchair inside a vehicle, this must have a securement system installed.

WARNING!

If, for any reason, it is not possible to transfer the user of the wheelchair onto a seat of the vehicle, the wheelchair can be used as a seat in a vehicle as long as the following rules and procedures are followed:

- The wheelchair must be fixed in the vehicle using a four-point securement system for wheelchairs.***
- The user must wear a three-point securement system for passengers fixed to the vehicle.***
- In the wheelchair the user must have the posture support belt fastened.***

Both the four-point securement system for wheelchairs and the three-point securement system for passengers have to be approved in compliance to the ISO-10542-2 regulation

WARNING!

Before setting off, contact the person responsible for transportation and request information on the availability of the equipment here listed.

Ensure there is enough free space around the wheelchair and the user to avoid the latter touching other passengers of the vehicle, non-padded parts of the vehicle, the accessories of the wheelchair and the anchoring points of the securement system.

WARNING!

Verify that the lifting points of the wheelchair are not damaged and that the parking brakes work perfectly.

The use of anti puncture tyres is recommended for transportation to avoid problems with the brakes due to the reduction in pressure of the tyres.

WARNING!

Injuries or damages could verify themselves due to components or accessories of the wheelchair loosening after an impact or an abrupt stop. Ensure that all removable or extractable components and accessories are removed from the wheelchair and stored safely in the vehicle. It is fundamental that the wheelchair is checked by a specialist after an accident, collision etc.

IMPORTANT!

It is recommended to provide training on the correct use of a product with a securement system. Refer to the user's manual supplied with the securement system. The following pictures can differ depending on the securement system's supplier.

Wheelchair fixing by means of a 4-point securement system.

CAUTION!

Inside the vehicle, place the wheelchair occupied by its user facing towards the direction of travel. Activate the wheelchair parking brakes. Activate the anti-tip system (if installed).

The anchoring positions of the wheelchair, where the securement system belts have to be placed, are tagged with symbols showing snap-hooks.

1. When using both front and rear belts of the 4-point securement system, fix the wheelchair to the snap-hooks assembled in the vehicle. Refer to the User Manual provided with the 4-point securement system.

Anchoring positions for the front belts

1. Fix the front belts onto the wheel stands.
2. Fasten the front belts to the snap-hooks system according to the instructions suggested by the safety belt supplier.
3. Release the parking brakes and fix the front tensioning belts pulling the wheelchair backwards.
4. Reactivate the parking brakes.

Anchoring positions for the snap-hooks at the rear

1. Fix the snap-hooks to the carbon fibre axle.
2. Fix the rear belts to the sideboard system according to the instructions suggested by the securement belt supplier.
3. Tighten the belts.

IMPORTANT!

Make sure the snap-hooks are covered by non-slip material in order to avoid side sliding of the axle.

IMPORTANT!

Check that the pivots are completely inserted on both sides and placed in the same position as the sideboard section. Make sure that inclination between the snap-hooks and the belts is between 40° and 45°.

Adjustment of the posture support belt

CAUTION!

The posture support belt must be used in addition, but it can never replace a homologated safety system for passengers (3-point belt).

1. Adjust the posture maintenance belt so that it adapts to the wheelchair user.

CAUTION!

Make sure the 3-point securement system for passengers is fixed to the user's body as tight as possible, without causing constraint to the user nor twisting of the belt.

Make sure the 3-point securement system is not kept away from the user's body by any part of the wheelchair such as armrests or wheels, etc.

Make sure the user's safety system goes from the user to the anchoring point without any interference by any other part of the vehicle, wheelchair, seats or accessories.

Make sure the sub-abdominal belt passes exactly on top of the user's pelvis and that it does not slide up towards the abdominal zone.

Make sure the user is able to reach the release mechanism without any need for help.

IMPORTANT!

Apply the pelvic part of the 3-point securement system well along the pelvis so that the inclination of the pelvic belt curves into the preferred zone (A) between 30° and 75° respect to the horizontal zone. A higher inclination is preferred, but never above 75°.

8.MAINTENANCE

8.1 SAFETY INDICATIONS

CAUTION!

Some materials are subject to natural ageing. That could cause some parts of the wheelchair to be damaged.

Get your wheelchair checked by a specialised dealer after a prolonged period or at least once a year.

Your Aria Wheels Srl dealer will be glad to assist you in carrying out regular maintenance interventions. In order to find an authorized dealer in the proximity, please get in touch with the Aria Wheels Srl contact at the end of the Manual.

8.2 MAINTENANCE PLAN

In order to guarantee the needed safety and reliability, carry out regularly the following maintenance interventions or have another person to carry them out.

	Weekly	Monthly	Yearly
<i>Check tyre pressure</i>	X		
<i>Check correct position of rear wheels</i>	X		
<i>Check adjustment of the backrest joint</i>	X		
<i>Check the posture maintenance belt</i>	X		
<i>Visual check</i>		X	
<i>Cleaning of front wheels</i>		X	
<i>Check bolt connections</i>		X	
<i>Check wheel spokes</i>		X	
<i>Check parking brakes</i>		X	
<i>Get wheelchair checked by a specialized dealer</i>			X

Check tyre pressure

Measure the tyre pressure checking with the relative pressures, (see the chapter 11.3 Tyres).

1. Inflate the tyre to the required pressure.
2. Check the tyre profile as well.
3. If necessary, replace the tyres.

Check the correct position of the rear wheels

1. Pull the rear wheel to check that the axle is correctly positioned in its location. It should not be possible to extract the wheel.
2. In case the rear wheels do not block into their location correctly, remove any deposit of dirt and scale. If the problem persists, get a specialized dealer to re-adjust the axle.

Check the adjustment of the backrest joint

1. Sit in the wheelchair and lean back onto the backrest. The backrest must be correctly blocked.
2. Pull the release cord. The pivots must be able to move freely on both sides. Otherwise, the backrest joint will have to be re-assembled by a specialized dealer.

Check the posture maintenance belt

1. Check that the posture maintenance belt is adjusted correctly.

IMPORTANT!

***Loose posture maintenance belts must be adjusted by a specialized dealer.
Damaged posture maintenance belts must be replaced by a specialized dealer.***

Visual check

1. Check if the wheelchair presents any loose parts, cracks or other defects.
2. In case any defects are identified, get a specialized dealer to immediately check the wheelchair.

Cleaning the front wheels

1. Check if the front wheels can be rotated freely.
2. Remove any dirt and hair from the front wheel bearings.

Check on bolt tightening

Bolts can loosen due to a constant use of the wheelchair.

1. Check that bolts are correctly tightened (footrest, seat cover, sides, backrest, frame, seat module).
2. Tighten any loose nuts with an adequate torque.

IMPORTANT!

***Safety nuts and bolts lose their efficiency after they have been loosened and tightened several times.
Get a specialized dealer to replace the safety nuts and bolts.***

Check spoke tension

Spokes must not be loose nor deformed.

1. Get a specialized dealer to adjust the tension of loose spokes.
2. Get a specialized dealer to replace any broken spokes.

Check parking brakes

1. Check the correct position of the parking brakes in their location.
The brake is adjusted correctly if, when the brake is in use, the brake shoe penetrates into the tyre for a few millimetres.
2. In case you realise the adjustment is incorrect, get a specialized dealer to adjust the brake correctly.

IMPORTANT!

After replacing or changing the position of the rear wheels, it's necessary to re-adjust the parking brakes.

IMPORTANT!

*The wheelchair can suffer damage caused by a strong collision or impact and not visible to the naked eye.
It is therefore fundamental to get a specialist to check the wheelchair after a strong collision or impact.*

Repairing or replacement of the inner tube

1. Disassemble the rear wheel and empty the inner tube from any remaining air.
2. Lift a heel of the tyre from the edge of the wheel rim. Use a device to remove tyres from bicycles; don't use a sharpened object, for example a screwdriver, as to avoid damaging the inner tube.
3. Extract the inner tube from the tyre.
4. Repair the inner tube by using a repairing kit for bicycles or, if necessary, replace the inner tube.
5. Slightly inflate the inner tube until it achieves a rounded shape.
6. Introduce the valve into the valve hole in the wheel rim and insert the inner tube into the tyre (the inner tube should adapt to the circumference of the tyre without forming any creases).
7. Lift the heels of the tyre on the edge of the wheel rim. Start from the area around the valve and use a device to remove tyres from bicycles. Check all along the circumference that the inner tube is not stuck anywhere between the tyre and the wheel rim.
8. Inflate the tyre to its maximum pressure. Make sure there is no air leaking from the tyre.

8.3 CARE

Your wheelchair will accompany you for many years if you take care of it regularly.

IMPORTANT!

*Sand and sea water may damage the ball bearings and the steel parts may rust if the surface is damaged.
Expose the wheelchair to sand and sea water only for a short time and clean it after having been on the beach.*

IMPORTANT!

Don't use abrasive substances, aggressive cleansers and high pressure cleaning devices.

1. Clean the lining and the metal parts with a damp soft cloth.
2. Dry the wheelchair accurately after using it in the rain or in the shower.
3. If the wheelchair is dirty, soften the dirt and remove it; then dry the wheelchair with care.

8.4 DISINFECTION

Note the concentration and exposure times of the disinfectant manufacturer. The product is suitable for spray and wipe disinfection with ordinary household disinfectants. All surfaces should be wiped with a clean cloth moistened with disinfectant. Apply the disinfectant evenly. Do not rinse and allow the product to air dry. Subsequently, the product must be checked for cleanliness and damages.

IMPORTANT !

To receive a list of all homologated disinfecting products, please ask your trusted dealer.

9 BREAKDOWN SOLUTIONS

9.1 SAFETY INDICATIONS

Daily use, new regulations or variation of the forces the wheelchair is subjected to may cause breakdown. The table below provides indications on how to recognise and eliminate breakdowns.

CAUTION!

In case anomalies of the wheelchair are noticed, e.g. a change in speed behaviour, contact your specialized dealer immediately.

IMPORTANT!

Some of the reported interventions must be carried out by an authorized specialized dealer. Those have been duly marked. We suggest all adjustments are carried out by a specialized dealer.

BREAKDOWN	POSSIBLE CAUSE	SOLUTION
<i>The wheelchair does not move forward in a straight line</i>	<i>The pressure of one of the rear wheels is incorrect</i>	<i>Adjust the tyre pressure, >11.3 Tyres,</i>
	<i>One or more spokes are broken</i>	<i>Replace the broken spoke(s), > Specialized dealer</i>
	<i>Spoke tensions are different</i>	<i>Get the tension of the loose spokes adjusted, > Specialized dealer</i>
	<i>The front wheel bearing is dirty or damaged</i>	<i>Clean or replace the bearing, > Specialized dealer</i>
<i>The wheelchair tips backwards easily</i>	<i>Rear wheels have been assembled too forward</i>	<i>Adjust barycentre, > Specialized dealer</i>
	<i>Excessive inclination of the backrest</i>	<i>Reduce inclination of the backrest, > Specialized dealer</i>
<i>Brakes intervene incorrectly or asymmetrically</i>	<i>Tyre pressure of one or both rear wheels is incorrect</i>	<i>Adjust tyre pressure, >11.3 Tyres,</i>
	<i>Brake adjustment is incorrect</i>	<i>Correct brake adjustment, > Specialized dealer</i>
<i>Resistance to movement is very high</i>	<i>Tyre pressure of the rear wheels is insufficient</i>	<i>Adjust tyre pressure, >11.3 Tyres,</i>
	<i>Rear wheels aren't parallel</i>	<i>Check that rear wheels are parallel, > Specialized dealer</i>
<i>Front wheels swerve when movement speed is high</i>	<i>Traction of the front wheel bearing block is insufficient</i>	<i>Slightly block the bolt on the bearing block axle, > Specialized dealer</i>
	<i>The front wheel is smooth</i>	<i>Replace the front wheel, > Specialized dealer</i>
<i>Front wheel steers with difficulty or is blocked</i>	<i>The ball bearing is dirty or defected</i>	<i>Clean or replace the ball bearing, > Specialized dealer</i>

10. AFTER USE

10.1 REUSE

The wheelchair can be reused. It is necessary to carry out the following operations:

- Cleaning and disinfection according to descriptions in chapter 8.3 Care and chapter 8.4 Disinfection of this User Manual.
- Inspection according to description in chapter 8.2 Maintenance plan of this User Manual.
- Adjustment to the user, based on the Technical Assistance documentation available at Aria Wheels Srl main office.

10.2 DISPOSAL

Defend the environment providing an adequate disposal of your wheelchair. Observe the national law and local regulations for the disposal.

For a correct disposal, contact your specialized dealer or your council's public administration as to get the address of a local disposal centre.

11. TECHNICAL DATA

11.1 WEIGHT AND DIMENSIONS

All weight and dimension specifications refer to a width of the seat and to its depth in all the possible configurations of the wheelchair. Dimensions and weight may change according to the different configurations.

Seat width:	240 - 460 mm (in increment of 20 mm)
Seat depth:	280 - 460 mm (in increment of 20 mm)
Seat to footrest:	120 - 180 mm (in increment of 20 mm)
Backrest height:	240 - 420 mm (in increment of 15 mm)
Knee-to-heel length:	300 - 440 mm (in increment of 10 mm)
Front seat height:	450 - 520 mm (in increment of 10 mm)
Rear seat height:	360 - 430 mm (in increment of 10 mm)
Centre of gravity	100 - 150 mm (in increment of 10 mm)
Backrest angle:	85° - 95° to floor
Rear wheel axle:	0° or 3°
Frame angle	83° for Aria 1.0Mg, ULTRA and SPECIALE 83° o 93° available only for Aria 2.0Al
Wheelchair weight SW440 with standard rear wheels: (postural backrest included)	approximately 9.5 kg* for Aria 1.0 approximately 10.7 kg* for Aria 2.0 approximately 8.1 kg* for Aria ULTRA approximately 7.1 kg* for Aria SPECIALE
Wheelchair weight SW440 without rear wheels: (postural backrest included)	approximately 6.3 kg* for Aria 1.0 approximately 7.5 kg* for Aria 2.0 approximately 4.9 kg* for Aria ULTRA approximately 3.9 kg* for Aria SPECIALE
Nominal load:	120 kg Aria 1.0 and Aria 2.0 100 Kg ULTRA and SPECIALE 75 Kg for KID

The weight depends on the possible configurations. Some dimensions may not be available for some wheelchair models. See the product order forms for the specs.

11.2 ENVIRONMENTAL CONDITIONS

Do not expose the wheelchair to temperatures below -20 °C and above 40 °C.

11.3 TYRES

The optimal pressure depends on the kind of tyre:

Max Tyre Pressure

High ease of movement tyres 7 bar 700 kPa 101 psi

Profiled tyres (Marathon) 7,5 bar 750 kPa 108 psi

Slick tyres (speed run) 10 bar 1000 kPa 145 psi

Mountain bike tyres 4 bar 400 kPa 58 psi

Solid rubber tyres - -

TOTAL WIDTH (TW)
(camber 0°) SW + 180 mm
(camber 3°) SW + 225 mm

TOTAL DEPTH (TD)
83° = SD + 120mm + 130mm + 100mm + (100mm - CG)
93° = SD + 120mm + 130mm + 100mm + (100mm - CG) + 35mm

11.4 LABELS

The symbols mentioned on the label are those compliant with regulations.

SYMBOLS	MEANING
	Catalogue number
	Manufacturer
	Date of Manufacture
	Read User Manual
	Serial number
	Compliant with CE regulations
	Rated flow
	Warning

12. DECLARATION OF CONFORMITY

Visit ariawheels.com to get the declarations of conformity

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[illegible]

[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

[illegible]

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

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