Monoflex™ Chest Support

For users who need strong but comfortable abdominal or chest support, the MonoflexTM is the perfect dynamic solution. Available in two styles.

One-Piece Underarm Release

Single pad for more surface area and stretch, with front-pull adjustment.

Proprietary high-strength Swivel Buckle on both sides provides reliable, secure closure and two-sided adjustment.



Two-Piece Center Release

Reduces bulk under the arms and fits under lateral supports for added stability. Great for independent users. Center Swivel Buckle equalizes tension and allows pads to contour around the shape of the body.





High-strength elastic padding provides just the right amount of stretch for comfort, while stabilizing even those with strong spastic tone. Double-tapered shape distributes pressure across the midsection while maintaining free arm movement. When worn on upper abdomen, it can improve breath control and speech. Durable, machine-washable material designed for the most demanding users.

Part Numbers

Size	One-Piece Underarm Release	Two-Piece Center Release
XS	SH120XS	SH122XS
S	SH120S	SH122S
М	SH120M	SH122M
L	SH120L	SH122L
XL	SH120XL	SH122XL

End-Fitting Options

Choose end-fitting and add suffix -B1, -B2 or -B3* to the end of the part number.

Flat-Mount™ Cinch-Mount™ Grommet Strap







* One-piece Monoflex™ only

(IP) Bodypoint[®]