



**"WITHIN 1 HOUR OF DELIVERY I WAS IN THE BUSH AND ON THE BEACH! A DREAM COME TRUE AND PROBABLY ONE OF THE BEST FEELINGS I'VE EVER HAD, ABSOLUTE FREEDOM."**

**Nick**  
New South Wales

**magic**  
*Mobility*™

Extreme X8

**OFF-ROAD  
AND OFF  
THE MAP**

**HIT THE ROAD**

Or the sand, or the bush, or the snow.\*

\*Not to be driven through water

**MAKING YOUR CHAIR WORK FOR YOU**

Our wheelchairs are found on building sites and farms, in stables and national parks - allowing people to continue to work.



**DO YOU KNOW YOUR LIMITS?**

The Extreme X8 will help challenge you to embrace your freedom.



V6/AUST/NOV21

TAILORED FOR  
YOUR LIFESTYLE



AUSTRALIAN  
MADE

Extreme  
**X8**™





# MAGIC MOBILITY EXTREME X8

The Magic Mobility Extreme X8 provides the freedom you've always wanted. The X8 is your partner to almost any destination. Sand, snow and every terrain you can throw at it. With unrivalled off-road performance and reliability, the X8 is perfect for outdoor adventurers, explorers and rural workers. If the great outdoors is where you want to be, the Extreme X8 is waiting to take you there.

"WE'VE BEEN  
CREATING  
MADE-TO-ORDER  
WHEELCHAIRS  
FOR OVER 25  
YEARS, SO THAT  
YOU CAN LIVE THE  
LIFE YOU  
CHOOSE"

Ashley Daff  
R&D Director




## ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.

## QUALITY BUILT

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



## ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



## PROUDLY AUSTRALIAN DESIGNED AND MADE

We love that our product is made in Australia. We love meeting our customers, hearing your ever changing needs and taking it back to the drawing board.

### Basic specifications:\*

- Length: 1030mm
- Width: 710mm
- Seat to floor height: 470mm (490mm with swing away leg rests)
- Weight limit: 182kg (400lbs), 155kg (340lbs) with lift and/or tilt
- Top speed: 10kph (6.2mph)
- Motors: 700W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Gel Batteries
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Rehab seat pan (from 12"-24" wide).

### Rehab positioning:

- Power options: 300mm seat lift, tilt (45 degrees), power recline, lay down, power elevating legrests
- Seating and backrests: MPS, rehab backrest, adjustable seat pans and seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, side mount, swingaway, power or gas elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any specific requirement
- A full range of armpads, straps, pads, harness, restraints, supports are available.

### Options and accessories:

- One-click activator (CJSM 2 or JSM joystick)
- Fold forward back
- Puncture resistant tyre liners
- 24V LED Lighting kit
- 12V Accessory socket and 2 x USB charger
- Jack
- Luggage rack, bag hooks, bags, pockets, keep safe secret boxes, tables, oxygen and vent carriers
- Special controls: Attendant, head, chin, foot, sip and puff, environmental controls
- Crash test compliant tie downs
- RAM® mounting system (tablet and swing away tablet holder, mobile phone holder, swing away camera mount, attendant cup holder, action camera holder and cup holder).

### The right wheelchair for:\*\*

- Extreme off-road adventures
- Everyday off-road country living, working and rural productivity
- Active outer urban living.

### Base options:



FLAME RED - BLACK RIMS



OCEAN BLUE - SILVER RIMS



CAMO FOREST - BLACK RIMS

\*Specifications subject to change. Check online for current version  
\*\*Please remember that your wheelchair has electric motors and must not be driven through water, into rivers, creeks and the sea



Because the world's not flat



FAMILY OF BRANDS

Magic Mobility +61 3 8791 5600  
sales@magicmobility.com.au

Explore with us [magicmobility.com.au](http://magicmobility.com.au)



Scan for agent locations