

Aria Wheels Srls

Registered Office_Via della Divisione Torino, 92 00143 Rome - Italy
 Production Department_Via A.Volta 7/G 36030 Costabissara (VI) - Italy
 TEL. / FAX +39 04441322533 | VAT id IT 12848301003 | email : info@ariawheels.com



Ultralight adjustable wheelchair, the rigid frame is assembled with structural adhesives and it is made of magnesium alloy and Ergal components.

Standard features: Carbon fibre rigid postural backrest BKRTR + folding system, central axle in carbon fibre with Aria® patented shock absorbing technology, carbon fibre side guards, carbon fibre rigid seat pan, stainless steel screws + the STD optional.

Maximum user weight: 120 kg
 Transport weight ± 5.4 kg
 (see the symbols for the related configuration)



€ 5.600,00

WITH STANDARD FEATURES
 (VAT and shipping cost not included)

In some countries the price may change by 15%

CUSTOMER NAME AND SURNAME: mako QUOTATION ORDER

DEALER NAME: IF co.,ltd. DELIVERY ADDRESS: 2-11-14,Nishi23jyo-kita,Obihiro,Hok

PRODUCT SPECIALIST: Noritaka Naito ORDER NUMBER AND DATE: _____


Before filling the boxes below regarding the sizes , please see the technical scheme in the last page.
Aria Wheels Srls in order to avoid mistakes on the order form, suggests to contact a qualified person for taking the measurements. Aria Wheels can not make changes to the product after receiving the order form.
 If no choices are made in the order form, the product will be sent in the standard configuration.


CUSTOMIZATION


Select your preferences with a flag in the box

SEAT WIDTH (SW) <i>Considered external tube to external tube of the frame</i>			NOT ADJUSTABLE AFTER SALE
<input checked="" type="radio"/>	SW340	Seat width 340 mm 13-1/4 inches	<p>Between the internal faces of the side guards: +10mm</p>
<input type="radio"/>	SW360	Seat width 360 mm 14-3/4 inches	
<input type="radio"/>	SW380	Seat width 380 mm 15 inches	
<input type="radio"/>	SW400	Seat width 400 mm 16 inches	
<input type="radio"/>	SW420	Seat width 420 mm 17 inches	
<input type="radio"/>	SW440	Seat width 440 mm 17-1/4 inches	
<input type="radio"/>	SW460	Seat width 460 mm 18 inches	

SEAT DEPTH (SD) <i>Considered as the depth of the carbon fibre seat plate</i>			NOT ADJUSTABLE AFTER SALE
<input type="radio"/>	SD360	Seat depth 360 mm 14-3/4 inches	
<input type="radio"/>	SD380	Seat depth 380 mm 15 inches	
<input checked="" type="radio"/>	SD400	Seat depth 400 mm 16 inches	
<input type="radio"/>	SD420	Seat depth 420 mm 17 inches	
<input type="radio"/>	SD440	Seat depth 440 mm 17-1/4 inches	
<input type="radio"/>	SD460	Seat depth 460 mm 18 inches	

SEAT TO FOOTREST (StF)				Considered from the front end of the seat plate to the tube of the footrest	NOT ADJUSTABLE AFTER SALE
<input type="radio"/>	StF120	Seat to footrest 120mm 4-3/4 inches	STD	※フロントフレーム延長 	◇
<input type="radio"/>	StF140	Seat to footrest 140mm 5-1/2 inches			
<input checked="" type="radio"/>	StF160	Seat to footrest 160mm 6-1/2 inches			
<input type="radio"/>	StF180	Seat to footrest 180mm 7 inches			


CAMBER REAR WHEELS			※キャンバー角度
<input checked="" type="radio"/>	CBR0	Camber 0°	
<input type="radio"/>	CBR3	Camber 3°	

CENTER OF GRAVITY (CG)				Considered from the backrest foam to the axle of the rear wheels	納車後調整可能
<input checked="" type="radio"/>	CG100	Center of gravity 100 mm 4 inches	if has been selected SD360 and StF120 this is the only one possible		※車軸前出し位置 
<input type="radio"/>	CG110	Center of gravity 110 mm 4-1/4 inches			
<input type="radio"/>	CG120	Center of gravity 120 mm 4-3/4 inches			
<input type="radio"/>	CG130	Center of gravity 130 mm 5-1/4 inches			
<input type="radio"/>	CG140	Center of gravity 140 mm 5-1/2 inches			
<input type="radio"/>	CG150	Centro di gravità 150 mm 6 inches			


SEAT HEIGHT REAR (SHR)				Considered from the floor to the back of seat plate	納車後、一定範囲内で調整可能
<input type="radio"/>	SHR360	Seat height rear 360 mm 14-1/4 inches	adjustments range -lower red frame- adjustments range -standard frame-	N.b.: 調整可能な範囲 the minimum dimensions are to be considered increased by ± 12mm SHR390 and 360 will be only possible with CG ≥ 120 ※後座高 別紙に注釈有り	◇
<input type="radio"/>	SHR370	Seat height rear 370 mm 14-1/2 inches			
<input type="radio"/>	SHR380	Seat height rear 380 mm 15 inches			
<input type="radio"/>	SHR390	Seat height rear 390 mm 15-1/4 inches			
<input checked="" type="radio"/>	SHR400	Seat height rear 400 mm 15-3/4 inches			
<input type="radio"/>	SHR410	Seat height rear 410 mm 16-1/4 inches			
<input type="radio"/>	SHR420	Seat height rear 420 mm 16-1/2 inches			
<input type="radio"/>	SHR430	Seat height rear 430 mm 17 inches			


SEAT HEIGHT FRONT (SHF)				Considered from the floor to the front of seat plate	納車後、一定範囲内で調整可能
<input checked="" type="radio"/>	SHF450	Seat height front 450 mm 17-3/4 inches	adjustments range -lower red frame- adjustments range -standard frame-	N.b.: 調整可能な範囲 the minimum dimensions are to be considered increased by ± 12mm SHR390 and 360 will be only possible with CG ≥ 120 ※前座高 別紙に注釈有り	◇
<input type="radio"/>	SHF460	Seat height front 460 mm 18-1/4 inches			
<input type="radio"/>	SHF470	Seat height front 470 mm 18-1/2 inches			
<input type="radio"/>	SHF480	Seat height front 480 mm 19 inches			
<input type="radio"/>	SHF490	Seat height front 490 mm 19-1/4 inches			
<input type="radio"/>	SHF500	Seat height front 500 mm 19-3/4 inches			
<input type="radio"/>	SHF510	Seat height front 510 mm 20 inches			
<input type="radio"/>	SHF520	Seat height front 520 mm 20-1/4 inches			

DURING THE ASSEMBLING SHR AND SHF DETERMINE IF WILL BE USED A STANDARD OR LOWERED FRAME

KNEE TO HEEL LENGTH (KHL) <i>Considered from the footrest cover to the the seat pan</i>			納車時に調整します
<input type="radio"/>	KHL300	Knee to heel length 300 mm 11-3/4 inches	<p>In the case of Flip-Back footrest the maximum KHL measurements are reduced of 40mm.</p>  <p>※フット長</p> <p>With 5" front wheels KHL400 needs SHF>460 KHL410 needs SHF>470</p> <p>With 5" front wheels KHL430 needs SHF>490 KHL440 needs SHF>500</p>
<input type="radio"/>	KHL310	Knee to heel length 310 mm 12-1/4 inches	
<input type="radio"/>	KHL320	Knee to heel length 320 mm 12-1/2 inches	
<input type="radio"/>	KHL330	Knee to heel length 330 mm 13 inches	
<input type="radio"/>	KHL340	Knee to heel length 340 mm 13-1/4 inches	
<input type="radio"/>	KHL350	Knee to heel length 350 mm 13-3/4 inches	
<input type="radio"/>	KHL360	Knee to heel length 360 mm 14-1/4 inches	
<input checked="" type="radio"/>	KHL370	Knee to heel length 370 mm 14-1/2 inches	
<input type="radio"/>	KHL380	Knee to heel length 380 mm 15 inches	
<input type="radio"/>	KHL390	Knee to heel length 390 mm 15-1/4 inches	
<input type="radio"/>	KHL400	Knee to heel length 400 mm 15-3/4 inches	
<input type="radio"/>	KHL410	Knee to heel length 410 mm 16-1/4 inches	
<input type="radio"/>	KHL420	Knee to heel length 420 mm 16-1/2 inches	
<input type="radio"/>	KHL430	Knee to heel length 430 mm 17 inches	
<input type="radio"/>	KHL440	Knee to heel length 440 mm 17-1/4 inches	

TYPE OF CARBON FIBRE SEAT PAN ※座シート			動きがある仕様とする	NEW
<input type="radio"/>	CARB2.3	Rigid carbon fibre seat pan	STD	
<input checked="" type="radio"/>	CARB2.5	Carbon fibre seat pan with transversal cuts for greater flexion		230,00 €

BACKREST HEIGHT (BRH) <i>Considered from the seat pan to the top of the backrest</i>			※背高さとバックレスト選択
<input type="radio"/>	BRH240	Backrest height 240 mm 9-1/2 inches	<p>調整可能な範囲</p> <p>※背高さと、装着されるバックレストサイズ</p>  <p>S様仕様はSize3で高さは390mmです</p> <p>for the dimensional specifications of the ARIA® postural backrests click here</p>
<input type="radio"/>	BRH255	Backrest height 255 mm 10 inches	
<input type="radio"/>	BRH270	Backrest height 270 mm 10-3/4 inches	
<input type="radio"/>	BRH285	Backrest height 285 mm 11-1/4 inches	
<input type="radio"/>	BRH300	Backrest height 300 mm 11-3/4 inches	
<input type="radio"/>	BRH315	Backrest height 315 mm 12-1/2 inches	
<input type="radio"/>	BRH330	Backrest height 330 mm 13 inches	
<input type="radio"/>	BRH345	Backrest height 345 mm 13-1/2 inches	
<input checked="" type="radio"/>	BRH360	Backrest height 360 mm 14-1/4 inches	
<input type="radio"/>	BRH375	Backrest height 375 mm 14-3/4 inches	
<input type="radio"/>	BRH390	Backrest height 390 mm 15-1/4 inches	
<input type="radio"/>	BRH405	Backrest height 405 mm 16 inches	
<input type="radio"/>	BRH420	Backrest height 420 mm 16-1/2 inches	

BACKREST ANGLE (BRA) <i>Considered as the angle than the floor</i>	※背パイプ角度
Indicate the desired angle between 85° and 95°	
90°	

FOOTREST <i>All models allow the KHL adjustment</i> ※フットレスト				NEW
<input checked="" type="radio"/>	FOOTA1	Inclinable carbon fibre plate	STD	
<input type="radio"/>	FOOTB1	Flip-Back with carbon fibre cover		230,00 €
<input type="radio"/>	FOOTF1	Fixed with carbon fibre cover		215,00 €

BRAKES ※ブレーキ				<2020/8/28変更> 操作性を重視
<input type="radio"/>	BRKASL	Aria® Scissor lock <small>not available if SD360 and StF120 have been selected</small>	STD	
<input checked="" type="radio"/>	BRKOPL	Out-Front Push to lock		110,00 €
<input type="radio"/>	BRKOSL	Out-Front Scissor lock		150,00 €
<input type="radio"/>	BRKAEL	Aluminium Eagle brakes		190,00 €

BRAKES SUPPORT ※ブレーキ取付ステー				NEW
<input checked="" type="radio"/>	BRKASL2	Standard clamp	STD	
<input type="radio"/>	CARB.U3	Carbon fibre clamp <small>not available in combination with Out-Front Push to lock code BRKOPL</small>		175,00 €

FRAME COLOR <i>Powder coating</i> ※フレームカラー							
<input type="radio"/>	FRMCW	White	STD	<input type="radio"/>	FRMCBL	Blue	180,00 €
<input type="radio"/>	FRMCG	Graphite grey	STD	<input type="radio"/>	FRMCY	Yellow	180,00 €
<input checked="" type="radio"/>	FRMCB	Black	STD	<input type="radio"/>	FRMCFR	Red	180,00 €
<input type="radio"/>	FRMGM	Metallic grey	STD	<input type="radio"/>	FRMCC	RAL _____	360,00 €

REAR WHEELS ※リヤホイール							
<input type="radio"/>	Without rear wheels <small>the width TW could be different from the one indicated on the last page</small>			- 200,00 €			
<input checked="" type="radio"/>	24" ↓	<input checked="" type="radio"/>	25" ↓				
<input type="radio"/>	RWHAG	<input type="radio"/>	RWHAG.1	Black Wheel with black hub and black spokes	STD		
<input type="radio"/>	RWHSG	<input type="radio"/>	RWHSG.1	Golz Six-Star with black hub and black spokes		675,00 €	
<input checked="" type="radio"/>	RWHTS	<input type="radio"/>	RWHTS.1	Golz Twin-Star Exchange with black hub		1.000,00 €	
特注	SELECT THE COLOR OF SPOKES →		<input type="radio"/> Black spokes	70,00 €	<input type="radio"/> Blue spokes	70,00 €	
			<input type="radio"/> White spokes	70,00 €	<input checked="" type="radio"/> Orange spokes	70,00 €	
<input checked="" type="radio"/>	RWHAL	<input type="radio"/>	RWHAL.1	Spinergy® Lite Extreme LX with black hub		810,00 €	
SELECT THE COLOR OF SPOKES →		<input type="radio"/> Black spokes	STD	<input type="radio"/> Red spokes	70 €	<input type="radio"/> Blue spokes	70 €
		<input type="radio"/> White spokes	STD	<input type="radio"/> Yellow spokes	70 €	<input type="radio"/> Orange spokes	70 €
				<input type="radio"/> Green spokes	70 €	<input type="radio"/> Purple spokes	70 €

HANDRIMS * <i>Not available with rear wheels Golz Six Star (RWHSG / RWHSG.1)</i> ※ハンドリム				
<input type="radio"/>	HDRMB	Black anodized aluminum	STD	
<input type="radio"/>	HDRMS	Silver anodized aluminum*	STD	
<input type="radio"/>	HDRMT	Aluminum with foam coated for tetraplegic*		130,00 €
<input type="radio"/>	HDRMTI	Titanium*		495,00 €
<input type="radio"/>	HDRFCFC	Carbo Life Curve L* <small>click here for details</small>		380,00 €
<input checked="" type="radio"/>	HDRFCFG	Carbo Life Gekko* <small>click here for details</small>		500,00 €
<input type="radio"/>	HDRFCFQ	Carbo Life Quadro* <small>click here for details</small>		500,00 €

REAR WHEELS TYRES ※タイヤ			
<input type="radio"/>	WTSWR	Schwalbe Right Run - Black/Grey	STD
<input type="radio"/>	WTSHX	Solid polyurethane	95,00 €
<input type="radio"/>	WTSWM	Schwalbe Marathon Plus Evolution - Black	95,00 €
<input checked="" type="radio"/>	WTSWO	Schwalbe One - Black	125,00 €

FORKS ※キャスターフォーク			
* With this selection, only 4" wheels code FRWH4 and 5" wheels code FRWH5 can be mounted.			
<input type="radio"/>	FRKDS2	Aria® Double - Sided - black anodized aluminum	STD
<input type="radio"/>	FRKMNA	Aria® Single - Sided - black anodized aluminum	195,00 €
<input type="radio"/>	FRKFC	Aria® Single - Sided - in carbon fibre	310,00 €
<input checked="" type="radio"/>	SMTH	Aria Smooth® - Shock absorbing forks* s//ooth®	465,00 €

2020/8/28決定♪

FRONT WHEELS ※キャスター車輪			
<input type="radio"/>	FRWH3	3" Black anodized aluminium wheels	STD
<input checked="" type="radio"/>	FRWH4	4" Black anodized aluminium wheels	STD
<input type="radio"/>	FRWC4	4" Carbon fibre wheels	180,00 €
<input type="radio"/>	FRWH5	5" Black anodized aluminium wheels	STD
5" front wheels can be selected only if the following parameters have been selected: 25" rear wheels SHR> 380 SHF> 470			

SIDE GUARDS			
<input type="radio"/>	CARB3	※固定式カーボンファイバー Carbon Fiber Rigid	STD
<input type="radio"/>	CARB3.1	Carbon Fiber Fendered ※フェンダー（泥除け）付き	290,00 €

SIDE GUARDS OPTION			
<input type="radio"/>	SDGPO	pull-out	120,00 €

※着脱式サイドガード

ARMRESTS ※アームレスト			
<input checked="" type="radio"/>	ARMRT	Non-removable armrests with PU coating	220,00 €
The armrests can be ordered exclusively in combination with the rigid side guards CARB3			

ARIA® POSTURAL BACKRESTS ※Aria バックレスト			
<input type="radio"/>		Without postural backrest and with Aria® clamp mounted code ALCMT15	- 200,00 €
<input type="radio"/>	BKRTR	Aria® rigid Postural Backrest - carbon fibre	STD
<input checked="" type="radio"/>	BKRTS	Aria® semi-rigid moisture-wicking Postural Backrest - carbon f.	410,00 € 320,00 €


see the dedicated order form by clicking here



BKRTR



BKRTS



integrated into the frame
ALCMT15

SEAT CUSHION ※シートクッション			
<input type="radio"/>	IMBOTT7	Cushion made of foam with breathable textile	85,00 €

PUSH HANDLES ※介助ハンドル (手押しハンドル)			
<input checked="" type="checkbox"/>	PHY	Adjustable push handles - black matte painted aluminum	110,00 €

BACKPACK ※背シート裏に背負うバッグ		ATTACHMENT ※スマートドライブ用アタッチメント	
<input type="radio"/>	BKPY	Cordura and textile	130,00 €
<input checked="" type="checkbox"/>	SMDC	For Smart drive	185,00 €

UNDER SEAT BAG / HANDLE			
<input checked="" type="radio"/>	USBG	Black cordura under seat bag with velcro fabric ※シート下のバッグ	60,00 €
<input type="radio"/>	MSSU	Under seat grab handle ※シート下の取って追加	65,00 €

ANTI TIPPER ※転倒防止装置			
<input type="radio"/>	ATPR	Anti tipper mounted on left	145,00 €
<input type="radio"/>	ATPR	Anti tipper mounted on right	145,00 €
<input type="radio"/>	ANTW	with wheels	45,00 €

CALF BAND ※足ベルト			
<input checked="" type="radio"/>	CLFB8	Calf band, height 8 cm - black textile	STD
<input type="radio"/>	CLFB6L	Calf band, height 6 cm - black leather	60,00 €




SCRATCH-RESISTANT PROTECTION ※フレームプロテクション (傷防止カバー)			
<input type="radio"/>	FRMP	Pair of frame protection made of black leather - velcro fabric	70,00 €

NOTES:


CODES - OPTIONAL AND ACCESSORIES



スポーク組み換え特注仕様



別紙に概略寸法図を添付します



The technical scheme do not consider the dimensions of any postural and anti decubitus cushions.

SW: Considered external tube to external tube of the frame. The dimension between the internal faces of the side guards is + 10 mm than the SW.

BRH: Considered from the seat plate to the top of the backrest.

SD: Considered as the depth of the seat plate. The seat plate can be slided of 20 mm in order to follow the different angles of the backrest.

CG: Considered from the backrest foam to the axle of the rear wheels.

FD: Frame Depth = StF + SD - 10mm + 130mm

Backrest foam, thickness 30 mm.

For a right-sizing of the wheelchair is recommended:

① SW | Seat Width _ A space of ± 5 mm between the user basin and the side guards

② SD | Seat Depth _ A space of ± 50 mm between the user's popliteal fossa and the rigid carbon seat.



TOTAL WIDTH (TW)

With Camber 0° SW ± 175 mm

With Camber 3° SW ± 220 mm

If you have selected the pull-out option on the side guards, add 6mm on the total width



TOTAL DEPTH (TD)

STANDARD FRAME:

StF + SD + 130 + 100 + (100-CG)

LOWERED FRAME:

StF + SD + 130 + 100 + (100-CG) - 25



TRANSPORT WEIGHT

Without rear wheels

± 5.4 kg

DO YOU WANT TO CHECK THE OVERALL SIZE OF YOUR CHAIR?
CLICK HERE AND FILL IN THE DATA, AUTOMATICALLY YOU WILL GET THE OVERALL MEASURES