



"WE HAVE A SIX ACRE PROPERTY AND THE EXTREME X8 ALLOWS ME TO GET TO ALL THE UNEVEN TERRAINS AND SEE PARTS OF THE FARM THAT I HAVE NOT BEEN ABLE TO ACCESS IN FOUR YEARS! THANK YOU FOR GIVING ME MY FREEDOM AND INDEPENDENCE."

**Michael Trost**  
Maryland USA

### HIT THE ROAD

Or the sand, or the bush, or the snow.\*

\*Not to be driven through water

### MAKING YOUR CHAIR WORK FOR YOU

Magic Mobility wheelchairs can be found on building sites, farms, stables, national parks - all allowing people to continue to work.



### DO YOU KNOW YOUR LIMITS?

The Extreme X8 will help you find them, challenge yourself and embrace your freedom.



# MAGIC MOBILITY EXTREME X8

AUSTRALIAN  
MADE  
-  
BUILT TO LAST  
-  
CUSTOMIZED  
FOR YOUR  
LIFESTYLE



201 Growth Parkway  
Angola, IN 46703  
United States

Tel 1-260-665-2769

Email [iim@mobility-usa.com](mailto:iim@mobility-usa.com)

Web [mobility-usa.com](http://mobility-usa.com)



# MAGIC MOBILITY EXTREME X8



The Extreme X8 is Magic Mobility's flagship powerchair. Unrivalled off-road performance and reliability make it perfect for outdoor adventurers, explorers and rural workers. The Extreme X8 is the most compact 4WD powerchair on the market.

#### Basic specifications:

- Length: 45½"
- Width: 28"
- Seat to floor height: 18½" (19½" with swingaway legrests)
- Weight limit: 400lbs 340lbs with lift and tilt
- Top speed: 6.2mph
- Motors: 700W Gear in Line
- Wheels: Off-road 14" (145/70-6) black
- Batteries: Group 24, 75Ah LA/AGM and Gel
- Penny and Giles controller with LCD display (R-Net)
- Centre post legrest with one piece footplate
- Sling backrest and rehab seat pan (from 12"-23" wide)\*.

#### Rehab positioning:

- Power options: 12" seat lift, tilt (45 degrees), power recline, lie down, power elevating legrests
- Seating and backrests: MPS, rehab backrest and adjustable seat pans and custom seating
- Headrests and mounts: MPS, rehab
- Legrests: centre post, side mount, swingaway, power elevating
- Footplates: 1 or 2 piece fixed/flip up/angle adjustable, any custom requirement
- A full range of armpads, straps, pads, harness, restraints, supports are all available from our extensive catalogue.

#### Options and accessories:

- 24V LED Lighting kit
  - USB charger
  - Jack
  - Luggage rack, bag hooks, gun mounts, bags, pockets, keep safe secret boxes, tables, drinks holder, oxy and vent carriers
  - Special controls: Attendant, head, chin, foot, sip and puff, environmental controls all available
  - Crash test compliant tie downs.
- The right wheelchair for:**
- Extreme off-road adventures
  - Everyday off-road country living, working and rural productivity
  - Active outer urban living.



SKY BLUE



BLACK



TANGERINE


"WE'VE BEEN  
CREATING CUSTOM  
MADE VEHICLES  
FOR 22 YEARS,  
SO THAT YOU CAN  
LIVE THE LIFE  
YOU CHOOSE"

**Ashley Daff**  
Director



#### ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know the small touches make a big difference.



#### ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in house engineering team are hands on, working side by side with customers and therapists. The entire Magic team work continuously on how we can improve the experience from choosing a Magic Mobility wheelchair to its performance in the places you want to go.



#### BUILT TO LAST

We never compromise on quality. With your wheelchair as your mobility, reliability is paramount. We'd much rather spend time improving our products than repairing them.



#### PROUDLY AUSTRALIAN OWNED AND DESIGNED. ASSEMBLED IN THE USA

We love that our product is made in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

\*Specifications subject to change.  
Check online for current version

\*Various model options excluded

**Magic**  
MOBILITY  
BECAUSE THE WORLD'S NOT FLAT

Visit us online to see how we change lives every day:  
[magicmobility.com.au](http://magicmobility.com.au)

**Postal Address**  
Magic Mobility Pty Ltd  
2/112 Browns Road  
Noble Park VIC 3174 Australia

**Tel** +61 3 8791 5600  
**Fax** +61 3 9701 0144  
**Email** [sales@magicmobility.com.au](mailto:sales@magicmobility.com.au)



Join the Magic Mobility social conversations!

